



### Features :

- Built with durable steel framework and high-quality ABS & PP plastic
- Three functions easily adjusted by remote control (back-rest, leg-rest and hi-low)
- Brake pedal is conveniently located below the foot board

### Technical Specifications:

Models	ToronCare 1030™	
	Dimensions	
Overall dimensions (LxWxH)	2160mmx1150mm x 445-655mm 85inx45.3in x 17.5 -25.8in	
Hi-Low	655mm-445mm 25.8-17.5in	
Bed Board dimensions (LxW)	1935mmx900mm 76.2inx35.4in	
Diameter of central brake castors	150mm 5.9in	
	Bed Articulations	
Back -rest	0-70°	
Leg -rest	0-35°	
	Weight	
Net Weight	112 Kg 247 lbs	
Load Capacity	<250 Kg <551 lbs	
	Accessories	
Standard Accessories (Included)	I.V Pole	