

Hydraulic Stretcher – ToronCare™ 1180



Features and Applications:

- Manual type transfer stretcher with double hydraulic system.
- 5 functions available by hydraulic, gas spring and crank system together.
- A pair of push-handles at the head and foot side, which help user to move the stretcher.
- Central locking castors are matched with direction castor under the chassis, helps the user for easy control of stretcher direction
- Collapsible stainless steel protective railing.
- Bed board for X-ray examination.



Technical Specifications:

Dimensions	
Overall dimensions (LxW)	2050x750mm (80.7x29.5 in)
Height adjustment (Bed board to ground)	850-625mm (33.5-24.6 in)
Size of bed board (LxW)	1930x630mm (75.9x24.8 in)
Diameter of Central-lock castor	150mm (5.9 in)
Articulations	
Back-rest lifting angle	0-60°
Leg-rest lifting angle	0-40°
Trendelenburg & Rev. Trendelenburg angle	0-13°
Weight	
Net Weight	105 Kg (231 lbs)



Standard Accessories (Included):

- I.V Pole

