

## Features and Applications:

- Manual type transfer stretcher.
- Single crank is used for height adjustment, stretcher can be positioned up and down by turning the crank clockwise and anti-clockwise respectively.
- Bright arrow symbols direct the user how to operate clearly.
- Central locking castors are matched with direction castor under the chassis, which helps the user for easy control of stretcher direction
- Patient on the mattress can be moved parallel to another bed or flat table by rotating the upper crank handle.
- Rotary rollover protective railing
- Insulating bed board



## Technical Specifications:

Dimensions	
<b>Overall dimensions (LxW)</b>	1955x635mm (76.9x25.0 in)
<b>Height adjustment (Bed board to ground)</b>	850-590mm (33.5-23.2 in)
<b>Size of mattress (LxW)</b>	1880x555mm (74.0x21.8 in)
<b>Diameter of Central-lock castor</b>	150mm (5.9 in)
<b>Horizontal move range</b>	400mm (15.7 in)
Weight	
<b>Net Weight</b>	92 Kg (203 lbs)



## Standard Accessories (Included):

- I.V Pole