

Features and Applications:

- Manual connecting Transfer Stretcher, can be used to defend cross-infection.
- Cart A is inside the operation room while Cart B is outside the operating room, both carts are connected by adjusting them to the same height and then pushing Cart B to Cart A, following the orbit part C which can be pushed from one cart to another and vice versa.
- Drop-down protective railing, PP bed produced by blowing molding.



Technical Specifications:

Dimensions	
Cart A & B connected (LxW)	3880x840mm (152.7x33.0 in)
Height adjustment (Bed board C to ground)	930-680mm (36.6-26.7 in)
Size of bed board C (LxW)	1905x610mm (75.0x24.0 in)
Articulation	
Back-rest lifting angle	0-60° ± 5°
Weight	
Net Weight	139 Kg (306 lbs)



Standard Accessories (Included):

- Mattress
- I.V Pole