



Features:

- Detachable wooden head and foot board
- Collapsible aluminum protective side rail with user-friendly control switch
- Powder-coated steel bed board with ventilated holes for increased comfort
- Bed equipped with 4 pc stand-board instead of castors, providing increased safety support for patient
- Three functions easily adjusted by remote control (back-rest, leg-rest and hi-low)

Technical Specifications:

Dimensions	
Overall dimensions (LxWxH)	2080mmx1006mm x 345-610mm 81.9in x 39.6in x 13.6-24.0in
Hi-Low	610mm-345mm 24.0in-13.6in
Bed Board dimensions (LxW)	1960mmx850mm 77.2in x 33.5in
Bed Articulations	
Back-rest	0-70°
Leg-rest	0-35°
Weight	
Net Weight	100 Kg 220 lbs

Standard Accessories (Included)

- IV Pole

