



### Features:

- Detachable wooden head and foot board
- Collapsible aluminum protective side rail with user-friendly control switch
- Powder-coated steel bed board with ventilated holes for increased comfort
- Back-rest and leg-rest are easily adjusted by remote control

### Technical Specifications:

Dimensions	
Overall dimensions (LxWxH)	2100mmx1135mm x 380mm 82.7in x 44.7in x 15.0in
Fixed height	380mm 15.0in
Bed Board dimensions (LxW)	1960mmx950mm 77.1in x 37.4in
Bed Articulations	
Back-rest	0-75°
Leg-rest	0-45°
Weight	
Net Weight	80 Kg 176 lbs

### Standard Accessories (Included)

- IV Pole

